ENGLISH

Indication:

The dental diamond tip is a device for use in dental treatment that is connected to high-rotation instruments (high-rotation handpiece) and eventually at low rotation (micromotor coupled in counter-angle with adapter chuck). It is used in several clinical situations, the most common being the realization of cavitary preparations, in addition to being used in the wear of a wide variety of materials found in dentistry. Being dental materials: restorative materials (hard metals, amalgam, porcelain, resinous composites, glass ionomers, in addition to dental tissues, such as enamel).

3R Indústria E Comércio Ltda is not responsible when the product is used outside the intended use proposed.

Use applications by type of dental diamond tips (sterile and non-sterile):

Dental Diamond Tips (sterile and non-sterile) can be used in high-rotation instruments (high-rotation pen), always with water/air cooling to reduce the heat generated by friction, as well as to eliminate residues of dental tissue and / or materials.

The indication of use and function of the device is related to the design and size of the active tip that will be used according to the clinical situations. The granulometry (grain size) and/or type of grain are important specifications on the diamond tips that guarantee the interchangeability and compliance of these products with the requirements of dental practice.

Dental diamond tips have different granulations: extra coarse (181 µm), thick (151 µm), medium (91 - 126 µm), fine (46 µm) and extrafine (30 µm). Extra thick (black) and thick (green) tips are used to quickly reduce the thick layers of the tooth structure. They promote a faster reduction than the middle tips, which, in turn, are used in routine tasks involving enamel, dentin and cavitary preparations in general. Fine tips (red color) and extrafine tips (yellow color) are used in the finishing of the wall and on the margins of the cavitary preparation in the surface finish of direct or indirect restorations of composite resin or ceramics.

The dental diamond burs (sterile and non-sterile) has different shapes and sizes, and each one is coded (ISO code) for identification.

A	B	C	D	E
806	31 + 4	001	524	006
Material of Working part (Diamond)	Shank and overall length (FG / RA)	Form (round, cylinder, etc.)	Type (Toothing, coarse, etc.)	Nominal size Largest working part diameter in 1/10mm
ISO 6360-1	ISO 1797-1 / ISO 6360-1	ISO 6360-2	ISO 6360-3	ISO 2157
COD. ISO:	806.314.001.524.006			

Source: Standard ISO 6360-1

The following dental diamond tips (sterile and non-sterile) are examples of tips commonly used in cavitary preparations:

- Spherical:** Indicated for the initial access to the caries lesion. Other uses are recommended, according to the diameter of the device:
 - Small Diameter:** minimally invasive cavitary preparation and additional retentions in cavitary preparation for restoration in silver amalgam.
 - Average Diameter:** delimitation of orientation grooves in preparations for total crown, facets and indirect restorations in general.

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03/May/2023

Page 2 of 13

- **Large Diameter:** coronary opening, removal of the roof of the pulp chamber and bevel making in the cavo-superficial angle of extensive cavitory preparations in direct and indirect restorations.
- **Spherical with collar:** Indicated for use in cavitory preparations, assisting, for example, in the initial access to carious lesions, opening or retention in direct dentistry, making orientation grooves in indirect restorations and regularization of the cavo-superficial angle.
- **Conic:** Although the shape of these diamond tips is conical, they have variations that will condition the clinical application, described below:
 - **Inverted Conic:** Mainly used to determine mechanical retentions. They are useful in the preparation of self-retentive cavitory preparations in extensive restorations. In addition, the larger ones flatten and regularize the gingival wall, pulp and dihedral angles.
 - **Conical Torpedo:** Its tip crosses the point of proximal contact in cavity preparations inside and outside the access. These tips are used on proximal faces for indirect restorations of metal, porcelain or composite resin.
 - **Conical Flat End:** Indicated to make cavitory preparations with expulsive walls and defined internal angles. It is also indicated for making niches to support the removable partial denture.
 - **Conical of rounded end:** Indicated to make cavitory preparation for inlay and onlay of porcelain / composite resin with expulsive walls and rounded angles and preparation for onlay with protection of cusps finished in conical chamfer. It is also indicated for making orientation grooves and reduction of the vestibular face in prosthetic preparations of laminates.
- **Flat-ended cylindrical:** They are used to perform parallelism on the walls of the preparation and make the dihedral angles sharper. The thin tips are used in the wall finish of the cavitory preparation. Indicated for making and finishing in the form of a shoulder or marginal step in prosthetic preparations.
- **Wheels:** Indicated to create orientation grooves for the wear of the vestibular face.
- **Flame:** They help in the making of concave bevel, preparations on the lingual/palatine faces in anterior teeth, finishing in the cavo-superficial angle and can also be useful in wear for occlusal adjustment.

Instructions for use (sterile and non-sterile)

- Insert the dental diamond tips only into a high-rotation pen that meets the requirements of the applicable IEC 60601 series of standards and other standards relevant to the product. Apply the abrasive part to the surface to be prepared. Use dental diamond tips only for the recommended indications.
- Select the dental diamond tips according to the speed (we have speed table - RPM- for use of the tips).

Diameter ISO	Speed Recommended	Maximum Speed	Work Pressure
(1/10mm)	RPM (Rotations per minute)	RPM (Rotations per minute)	Weight (g)
008-018	250.000 – 300.000	450.000	20 – 50
021-031	100.000 – 120.000	160.000	20 – 50

This review cancels and replaces the previous revision.
3R Indústria e Comércio Ltda
Rua Ptolomeu, 290 – Socorro
São Paulo - SP - CEP.: 04762-040 - Brasil

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Cleaning procedures (sterile and non-sterile diamond tips):

- Clean and sterilize dental diamond tips before each use as per cleaning and sterilization item
- It is recommended that cleaning be done immediately after the end of the procedure, to avoid drying out contaminants, which can make cleaning difficult.
- After cleaning, check the diamond tips to ensure that any traces of contamination have been removed. Repeat cleaning if necessary.
- The cleaning of dental diamond tips can be done in two ways: Manual or Ultrasonic.

Manual cleaning procedure

If manual cleaning is the only option available, do the following.

- Use an area specifically reserved for cleaning;
- Wash the dental diamond tips in running water and keep them immersed in a bath of neutral pH cleaning solution (enzymatic detergent) following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agente (enzymatic detergent). Care should be taken to avoid spreading contaminants by spraying during the brushing process. Do not use wire brushes.
- Care should be taken to clean cracks and other hard-to-reach areas. Inspect visually to confirm debris removal. Repeat the cycle if necessary.
- Wash dental diamond tips thoroughly under running water until they are visibly clean.
- Dry the device with clean compressed air or paper towel.

Ultrasonic cleaning procedure

For ultrasound cleaning, it is recommended to proceed as follows:

- Prepare a neutral pH cleaning solution (enzymatic detergent) following the agent manufacturer's instructions for correct concentration, exposure time, temperature, and water quality.
- Place the dental diamond tips in the appropriate place for the ultrasonic equipment, following the manufacturer's instructions for use.
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Wash dental diamond tips thoroughly under running water until they are visibly clean.
- Dry the device with clean compressed air or paper towel.

Sterilization Requirements for Dental Diamond Tips (sterile and non-sterile)

- This item applies to dental diamond tips purchased in non-sterile form and to the reuse of sterile purchased products.
- It is recommended to sterilize the Dental Diamond Tips by saturated steam autoclave under pressure at 135°C for 10 minutes, or the sterilization parameters should be followed according to the relevant policies/regulations of the country where the device is being used.
- It is the full responsibility of the professional to keep the autoclave equipment updated to ensure the correct sterilization of the product.

LIMITATIONS OF USE: The recommended shelf life of the product is 5 uses (5 reprocessing), but it is of paramount importance that before each use is performed the verification, that is if the dental diamond tips are in perfect conditions of use, otherwise they can overheat the dental tissues and lead to pulp response triggering painful sensation in the patient, in addition to not complying with the necessary wear requirement.

Contraindications

- To avoid fracture and distortion of the instrument, do not apply excessive force when inserting the dental diamond tips at high speed;
- Respect the working speed (rpm Max) defined on the product label, otherwise it may generate undesirable heat;
- Do not use dental diamond tips as incorrect levers and/or angles;
- Do not apply excessive pressure to the dental Diamond Tips, as this may cause undesirable heat and/or cause failures in the Tips and risks to the patient;
- Chlorine or chloride cleaning agents as an active ingredient are corrosive to stainless steel and should not be used.
- Prolonged immersion of dental diamond tips in disinfectant solutions can result in corrosion and should therefore be avoided.
- Do not use cold sterilization methods (chemical process) for sterilization of dental Diamond Tips. These chemical agents can cause corrosion of the products.

Warnings and Precautions:

- Always check that the sterilization is within the shelf life on the product label. If it is not within the expiration date, proceed according to the items of Cleaning and sterilization.

This review cancels and replaces the previous revision.

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03/May/2023

Page 4 of 13

- Sterilization will only be guaranteed if the packaging is not damaged, in case of damage do not use the product before proceeding according to the items of Cleaning and sterilization
- In the reuse of the sterile acquired Dental Diamond Tip, proceed according to the cleaning and sterilization items.
- Clean and sterilize dental diamond tips before each use as per cleaning and sterilization items
- For user and patient safety, the dental diamond tip should be used only by qualified dental professionals.
- Before use, make sure that the dental diamond tips are in perfect condition for use, especially in relation to the diamond area and integrity of the rod
- Before use, make sure that the rod of the dental diamond tips is fully seated in the high-rotation pen;
- Never force the dental diamond tips into the high-rotation pen. In case of difficulty of fitting, check the turbine and the dental diamond tips.
- Respect the working speed (RPM Max) defined on the product label, otherwise it may generate undesirable heat.

Using (sterile and non-sterile):

- Move the dental diamond tips continuously to avoid localized heating and consequent breakage of the tips causing risks to the patient;
- Use proper cooling methods to ensure proper cooling of the parts used for treatment;
- Use aspiration techniques that do not compromise the cooling of the parts involved in the treatment;
- Eye protection should be used to protect against ejected particles.
- Used dental diamond tips should be considered contaminated instruments and appropriate handling precautions are required. Wear personal protective equipment (PPE), including gloves, masks, and safety goggles when handling, as with all other dental instruments

Post-use (sterile and non-sterile):

- The reuse of the product causes wear of the abrasive surface reducing its performance.
- In cases of return of dental diamond tips to 3R (manufacturer), make sure that they have been effectively cleaned, sterilized and properly packed for return.
- LOT Number: It is important to keep the batch number, identified on the product label, to be informed of any occurrences with the manufacturer.
- The material must be disposed of in accordance with local regulations.

Storage conditions (sterile and non-sterile))

- 1.Storage should be dry and clean at room temperature.
- 2.Keep the product stored in its original packaging until the time of use.

Additional Information

- 1.The reprocessing processes described above have been validated and are suitable and effective for the reuse of dental diamond drills. 3R is not responsible if the instructions for use are not followed properly, resulting in undesirable results in the efficacy and safety of the product/patient.
- 2.If necessary, 3R provides instructions for use in printed form, at no additional cost, upon request by email (administrativo@3r.ind.br) or telephone (11-5525-0590 / 0800-773-5525), and we will return within a maximum of 7 working days from the patient's request.

Always check the website: www.microdont.com.br, the most current version of the usage instruction corresponding to the product.

Symbols:

Symbol	Symbol Title	Description
	Expiration / Use-by date	Indicates the date after which the health product should not be used.
	Consult instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date when the health product was manufactured
	Manufacturer	Indicates the Manufacturer of the product for the health as defined in the Directive of the European Community 90/385/EEC, 93/42/EEC and 98/79/EC.
	Batch Code	Indicates the manufacturer's batch code so that the batch or departure can be identified.
	Catalog number	Indicates the manufacturer's reference number so that the product for the health can be identified.
	Sterilized using irradiation	Indicates that the product is sterilized
	Non-Sterile	Indicates that the product is not sterile
	Recyclable	Indicates that packaging can be recycled
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	Authorized representative in the European Community	Indicates the Authorized representative in the European Community

This review cancels and replaces the previous revision.

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