INTERLIG

PORTUGUÊS

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ENGLISH

Braided glass fiber impregnated with light-cured composite resin

INTRODUCTION

Glass fiber was developed as a solution to combine strength and lightness. Although limited, its use in Dentistry started in the 1960's. With further advances in silane chemistry and adhesion technology, the 80's marked a true glass fiber revolution in a variety of fields such as Periodontics, Prosthodontics and Operative Dentistry. The esthetic characteristics of glass fiber along with its most desirable mechanical properties allow immediate, highly conservative, durable and low cost restorative and preventive dental work.

DEFINITION

Interlig is a structure of intertwined glass fibers impregnated with light-cured composite resin.

COMPOSITION

- Glass fibers (weight) - $60 \pm 5\%$.
- Impregnated resin (weight) $40 \pm 5\%$: Bis-GMA, diurethane, barium glass, silicon dioxide, catalysts.

CHARACTERISTICS

- Fiber structure: intertwined scaffold.
- Optical: translucent.
- Thickness: 0.2 millimeter (mm).
- Width: 2 mm.
- Length: 8.5 centimeters (cm).

INDICATIONS

- 1.Splinting in Periodontics, Orthodontics and of avulsed or extracted teeth.
- 2.Reinforcement of resin retained fixed bridges (with natural or artificial teeth) and temporary fixed bridges made directly (in the dental office) or indirectly (in the lab).
- 3.Reinforcement and repair of acrylic resin bridges.
- 4.Fabrication of space retainers in pediatric dentistry.

PRECAUTIONS

- The components of INTERLIG® may be allergenic. Avoid contact with oral tissues, skin and eyes. Use gloves, mask and eye protection. If accidental contact occurs, wash with water. If necessary, seek medical attention.
- Keep it in dry, cool place (below 28°C / 82.4°F), protected from light, eugenol products and excessive heat

(furnaces, sterilizers).

- Put the remaining fiber in the package immediately after use to avoid excessive exposure to environmental light.

TECHNIQUE FOR USE

A. PERIODONTAL SPLINT

- a. Pre-op;
- b. Isolate with rubber dam (ideally) and place wedges interproximally to help controlling the flow of resin in these areas;
- c. With a pencil, delineate the area for cutting a groove;
Cutting a groove is an optional procedure which depends on the clinical situation.
- d. Make a 0.8 mm deep x 2.0 mm wide groove;
- e. Etch with 37% phosphoric acid for 30 seconds; wash and dry (areas of exposed dentin should not be desiccated);
- f. Apply the adhesive FUSION-DURALINK® system (same procedures used for direct composite resin restorations);
- g. Light-cure all segments of the splint site for 20 seconds each;
- h. Fill half of the groove with a hybrid or microhybrid composite resin (microfil resins should not be used in splinting);
- i. Cut a section of INTERLIG® slightly shorter than the length of the groove; position INTERLIG® in the groove and press it towards the composite resin;
- j. Gently pressing, light-cure each segment of the splint for 20 seconds; curing of non-pressed areas should be avoided by covering them with an aluminum foil;
- k. Apply a second layer of composite resin to cover the groove totally; the resin must also be placed slightly beyond the limits of the groove;
- l. Light-cure each segment for 40 seconds;
- m. Remove the excess of resin, check the occlusion, and adjust it if necessary; proceed with finishing and polishing;
- n. Apply a final layer of fluid resin (adhesive) to fill any irregularity or superficial bubbles and voids; light-cure all segments for 20 seconds each;
- o. Final result.

B. DIRECT FIXED BRIDGE

- a. Pre-op;

- b. Mark the site where the root will be cut after extraction;
- c. Area after tooth removal;
- d. Cut the root 1 to 2 mm apically to the mark considering soft tissue shrinkage during the healing process; shaping and polishing of the apical end of the tooth should produce a pontic-like, smooth and ovate design;
- e. Open the canal through the root, remove all pulpal tissue and fill the canal with composite resin (phosphoric acid and an adhesive system must be used);
- f. Isolate with rubber dam and cut a groove (0.8 mm deep x 2.0 mm wide) on the abutment teeth and on the extracted tooth;
- g. Fix the extracted tooth on the desired site with a small amount of composite resin placed facially; do not acid-etch now because this composite resin will be removed later; light-cure the resin;
- h. Acid-etch (37% phosphoric acid) for 30 seconds; wash and dry without desiccating any exposed dentin;
- i. Apply the adhesive FUSION-DURALINK® system, wait for 20 seconds, and remove excess of bonding agent with a light stream of air; light-cure each segment for 20 seconds;
- j. Apply a layer of a hybrid or microhybrid composite resin filling half of the groove;
- k. Measure and cut INTERLIG® with a scalpel or scissors in a slightly shorter size than the groove to facilitate insertion;
- l. Place INTERLIG® in the groove, press it, and light-cure each segment for 40 seconds; curing of non-pressed areas should be avoided by covering them with an aluminum foil;
- m. Apply a new layer of composite resin to cover INTERLIG® completely and fill the groove; the resin must be placed a little beyond the limits of the groove; light-cure for 40 seconds;
- n. Remove excess of resin, check the occlusion, and adjust it if necessary; proceed with finishing and polishing; apply a layer of fluid resin (adhesive) to fill any irregularity or superficial bubbles and voids; light-cure all segments for 20 seconds each;
- o. Final result.

C. INDIRECT FIXED BRIDGE (FABRICATED ON A WORKING CAST)

- a. In the mouth, prepare a groove (0.8 mm deep x 2.0 mm wide) on the abutment teeth and take the impression;
- b. On the working cast, paint the groove and its adjacent areas with a thin layer of a release agent;
- c. To be used as a pontic, select an acrylic resin denture tooth or make a tooth with porcelain or composite resin;
- d. Apply the adhesive FUSION-DURALINK® system and light-cure for 20 seconds;
- e. Insert a layer of light-cured hybrid or microhybrid composite resin to fill half of the groove; do not cure;

- f. Measure and cut INTERLIG® with a scalpel or scissors in a slightly shorter size than the groove to facilitate insertion; place INTERLIG® in the groove;
- g. Press it and light-cure each segment for 40 seconds; curing of non-pressed areas should be avoided by covering them with an aluminum foil;
- h. Apply a new layer of composite resin to cover INTERLIG® completely and fill the groove; light-cure for 40 seconds;
- i. With carbide burs or discs, trim the stone to release the bridge; the restoration should not be forced out of the working cast to avoid the risk of fracturing the connections; proceed with the final adjustments and polishing;
- j. Bridge ready; it must be cemented with an adhesive technique (adhesive system and a composite resin or resin cement).

D.REINFORCEMENT OF TEMPORARY ACRYLIC RESIN BRIDGES

- a. Cut a 0.8 mm deep x 2.0 mm wide lingual groove on the bridge;
- b. Apply the adhesive FUSION-DURALINK® system on the groove and light-cure for 20 seconds;
- c. Measure and cut INTERLIG® with a scalpel or scissors in a slightly shorter size than the groove to facilitate insertion;
- d. Insert a layer of light-cured hybrid or microhybrid composite resin to fill half of the groove; do not cure;
- e. Place INTERLIG® in the groove;
- f. Press it and light-cure each segment for 40 seconds; curing of non-pressed areas should be avoided by covering them with an aluminum foil;
- g. Apply a new layer of composite resin to cover INTERLIG® completely and fill the groove; light-cure each segment for 40 seconds;
- h. Proceed with the final adjustments and polishing;
- i. Final result.

WARNINGS

- INTERLIG® must always be covered by the composite resin; any exposed fiber should be gently removed with burs and the area corrected with a new layer of composite resin.
- Structures reinforced with INTERLIG® should not create occlusal interferences; always check carefully.
- Harmful habits (grinding, clenching, biting objects) should be corrected before indicating restorations reinforced with INTERLIG®.

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