**INSTRUCTIONS FOR USE**  
**APPLICATORS – BRUSH**  
**APLICADORES - BRUSH**

**CE**  
**MS: 80676920019**

IFU-Applicators-Brush- Rev. 00  
 05/May/2023  
 Page 1 of 4

**ENGLISH**

**Indication:**

Disposable Applicators are used to apply: identification solutions (caries, plaque, etc.), adhesive agents, sealants, eugenol, hemostatic solutions, cavity liners, acid and fluoride conditioners. It has a bendable rod and small amounts of substances can be applied.

**Instructions for use:**

- Remove an Applicator from the tube.
- Touch the tuft of fibers in the solution to be applied.
- Apply the solution to the place to be treated.
- The professional can choose a color for each type of application, according to the size of the region to be worked, making it easier to identify the product.

**Contraindications:**

- Non-sterile product. 3R recommends single use.
- Product for professional use only.
- Keep out of the reach of children.
- Non-sterilizable product, discard after use.
- Do not subject the device to high temperatures.
- During operation, the professional must use all personal protective equipment, including protective eyewear.
- Wear protective eyewear on the patient.
- LOT number: It is important to keep the lot number, identified on the product label, to be informed in eventual occurrences with the manufacturer.
- Disposal of material must follow relevant local legislation.

**Warnings and Precautions:**





- Keep the product in its original packaging until use. Inspect the packaging before use.
- The product must be stored in an airy, dry place, protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.

**Other Additional Information:**

- If necessary, 3R provides the instructions for use in printed format, at no additional cost, upon request by e-mail (administrativo@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will return within a maximum of 7 working days from the request.
- Always check on the site: [www.microdont.com.br](http://www.microdont.com.br), the most up-to-date version of the usage instruction corresponding to the product.

**Composition:** Polypropylene and polyamide fibers.

**Symbology:**

Symbol	Symbol Title	Description
	Expiration / Use-by date	Indicates the date after which the health product should not be used.
	Non-Sterile	Indicates that the product is not sterile
	of a single use	do not reuse
	Consult instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date when the health product was manufactured
	CE Mark	Indicates that the product has CE Marking and the identification of the Notified Body. (if applicable)
	Keep away from sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the Manufacturer of the product for the health as defined in the Directive of the European Community 90/385/EEC, 93/42/EEC and 98/79/EC.
	Batch Code	Indicates the manufacturer's batch code so that the batch or departure can be identified.
	Catalog number	Indicates the manufacturer's reference number so that the product for the health can be identified.
	Authorized Representative	Authorized Representative

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