INSTRUCTIONS FOR USE

DENTAL CARBIDE BUR – NON-STERILE

PUNTAS DE CARBURO DENTAL – NO ESTÉRIL

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ENGLISH

Indication:

The carbide dental drill is a cutting device for use in dental treatments that, coupled to a micromotor of high or low rotation, is used to act on a wide variety of tissues, such as: bone, enamel, dentin and also restorative materials such as hard metals and composites.

a) Carbide Burs Special Types

- **Zekrya:** Surgical use, odontosection, root separation, removal of root remains, the indication of its size varies according to the region to be reached.
- **Endo-Z:** endodontic use, preparation/smoothing of the pulp chamber during access surgery.
- **Transmetal:** High cutting power, can remove restorative materials such as amalgam, in addition to the removal of metallic prosthetic parts such as metal crowns, cast metal restoration, among others.
- **Lindemann and Lindemann Cross Cut:** Surgical use, indicated for cutting hard/bone tissues. The cross-section ensures a more efficient cut.
- **Surgical Drills:** Used for surgical procedures, bone cuts and dental cuts.
- **Finishing Drills - 12 blades:** Used in the finishing of restorations and rounding of cavo-superficial angles, the choice of the drill is preferably made by the professional executor according to the region of the oral cavity and dental element.
- **Finishing Drills - 30 blades:** Used in finer and more delicate finishes in restorations and rounding of the cavo-superficial angle, it can also be used for the removal of orthodontic resin remains. The choice of drills is preferably of the professional executor according to the region in the oral cavity and dental element.

b) Carbide drills FG

- **Spherical:** Used for cavitory access, eventually for opening cavities in order to perform removal of carious tissue, delimitation of additional retention and finishing line for prosthetic preparations.
- **Pear:** It is used to make additional retentions, rounding of angles, and eventually finishing lines and cingulate finish.
- **Inverted Conic:** Cavitory preparations where the walls must converge to the occlusal and create retentions.
- **Cylindrical:** Access to regions of greater difficulty, finishing of parallel surrounding walls, pulp wall planning.
- **Conic:** Preparation for indirect restorations, realization of grooves and channels and eventually additional retentions.

c) Carbide drills AC

These drills are used in micromotor of low rotation (counter-angle), the indications are the same as the FG, but with smaller rotation, ideal for the removal of carious tissue near the pulp where the work must be meticulous to avoid the exposure of the same. Despite the lower rotation, it is still necessary to use air/water cooling to avoid the heat generated by the friction of the drill.

Fundamentals of its operation and action:

Carbide drills are mounted on low-speed or high-speed pens (dental turbines) that have high rotations and low torque, so they are small diameter grinding wheels. The drive of these instruments can be carried out by means of a compressed air turbine or an electric motor. Compressed air transmission generally provides lower available torque value and higher rotation. For electric motor drive, the available torque is higher; But the rotation is smaller. The Carbide Drills are produced in Tungsten Carbide and Stainless Steel, with measurement, tolerance, tensile strength and hardness.

Sterile carbide drill

Do not use if the package is open or damaged. The sterile carbide drill does not require prior cleaning and sterilization.

If you choose to reuse, follow the cleaning and resterilization guidelines outlined below.

Non-sterile – Carbide Drill

The non-sterile carbide dental drill requires prior cleaning and sterilization.

Follow the cleaning and sterilization guidelines outlined below before use.

Instructions for use:

1. Insert the Carbide Bur FG rod into the high-speed turbine attachment using the "button" or drill bag.
2. Or in the case of state-of-the-art AC carbide, unlock the fitting of the low-speed micromotor, insert the rod into the socket and lock the rod into the low-speed micromotor.
3. Select the carbide drill according to the procedure performed and follow the speed according to the table below.
4. Check that the rod is well seated in the micromotor and that it works properly.
5. Take the micromotor to the desired location.
6. Activate the micromotor.
7. Apply the active drill to the desired location intermittently.

This review cancels and replaces the previous revision.

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The product can be reused if it is properly cleaned and sterilized.

LIMITATIONS OF USE: Only 5 reprocessing is recommended.

The cleaning of the carbide drill can be carried out in two ways: Manual or Ultrasonic. The professional who performs the cleaning and sterilization of the device should be using all protective PPE to avoid the risk of cross-contamination.

Manual cleaning procedure:

- Use an area specifically reserved for cleaning.
- Rinse the dental carbide drills in running water and keep them immersed in a bath of enzymatic detergent solution, following the manufacturer's instructions.
- After soaking and keeping it submerged, brush well away from the body with the neutral cleaning agent. Care should be taken to avoid spreading contaminants from the spray during brushing procedures. Do not use wire brushes.
- Special care should be taken to clean cracks and other hard-to-reach areas. Inspect visually to confirm debris removal, repeat the cycle if necessary.
- Wash the carbide drill with running water until it is visibly clean.
- Dry the appliance with clean compressed air or paper towels.

Ultrasonic cleaning procedure

- Prepare an enzymatic detergent solution by following the agent manufacturer's instructions for the correct concentration, exposure time, temperature and water quality.
- Place the carbide drill in the appropriate place of the ultrasonic equipment, following the instructions for use of the equipment manufacturer.
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Wash the carbide drill with running water until it is visibly clean.
- Dry the appliance with clean compressed air or paper towels.

Sterilization requirements with carbide drills (sterile and non-sterile)

- This item applies to Carbide Bur tips purchased in non-sterile form and to the reuse of sterile purchased products.
- It is recommended that Carbide Bur be sterilized by steam autoclave saturated under pressure at 135°C for 10 minutes, properly packed in its own autoclave packaging or the sterilization parameters must be followed in accordance with the relevant policies/regulations of the country where the device is being used.
- It is the full responsibility of the professional to keep the autoclave equipment updated to ensure the correct sterilization of the product.

Contraindications:

- To avoid fracture and deformation of the instrument, do not apply too much force when inserting the carbide drill into the pens (dental turbine).
- Respect the working speed (Max rpm), otherwise it may generate unwanted heat.
- Do not use the carbide drill as a lever and/or at incorrect angles.
- During use, do not apply excessive pressure to the carbide drill, as this may generate unwanted heat and/or cause carbide drill failure and risk to the patient.
- Prolonged immersion of carbide drills in disinfectant solutions can cause corrosion and should therefore be avoided.
- Do not use cold sterilization methods (chemical process) to sterilize the carbide dental drill. These chemical agents can cause corrosion of the products.
- Do not use sodium hypochlorite as a cleaning agent, as it may cause corrosion in the product.
- Do not use the device in patients allergic to any tungsten carbide component.

Warnings and Precautions:

- Before use, check that the Carbide Drill is in perfect condition for use, especially in relation to the integrity of the stem.
- Before use, make sure that the drill bit is well fitted to the low- or high-speed micromotor and that the coolant is working properly.
- Never force the carbide drill into the turbine. In case of difficulty of assembly, check the turbine and the carbide dental drill.
- Respect the working speed (Max rpm), otherwise it may generate unwanted heat.
- Use proper cooling methods to ensure proper cooling of the parts used in the treatment.

This review cancels and replaces the previous revision.

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Page 3 of 10

MS: 80676920039

- Use suction techniques that do not compromise the cooling of the parts involved in the treatment.
- Safety goggles should be worn to protect against ejected particles.
- Wear personal protective equipment, including gloves, masks, safety goggles, and long-sleeved lab coats when handling the carbide drill, as you would with all other dental instruments. Wear goggles on the patient as well.
- In cases of return of the carbide drill, make sure that it has been properly cleaned, sterilized and packed for return.
- BATCH NUMBER: It is important to keep the batch number, identified on the product label, to be informed of possible incidents with the manufacturer.
- The disposal of the material must follow the relevant local legislation.

Conservation and Storage:

- For sterile supplied products: sterilization will only be guaranteed if the packaging is not open or damaged.
- The product must be stored in an airy, dry and protected from light place to maintain optimal storage and transport conditions, as well as physical and chemical integrity.
- Keep the product in its original packaging until the time of use. Inspect the packaging before use.
- The carrier must be informed of the contents. The care that must be taken in the transport is demonstrated in the packaging, thus ensuring the protection of the instrument from shipment to delivery to the customer. Packaging should not be left to the action of time.

Other additional information:

- The sterilization processes described above have been validated and are suitable and effective for the use of carbide cutters. 3R is not responsible if the instructions for use are not followed correctly, generating results not intended for the efficacy and safety of the product/patient.
- If necessary, 3R provides instructions for use in printed format, at no additional cost, upon request by e-mail (administrativo@3r.ind.br) or by telephone: (+55 11-5525-0590 / +55 0800-773-5525), and we will return it within a maximum of 7 working days from the request.
- Always consult the website: www.microodont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition: Tungsten carbide and stainless steel

Symbology:

Symbol	Symbol Title	Description
	Expiration / Use-by date	Indicates the date after which the health product should not be used.
	Consult instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date when the health product was manufactured
	Manufacturer	Indicates the Manufacturer of the product for the health as defined in the Directive of the European Community 90/385/EEC, 93/42/EEC and 98/79/EC.
	Batch Code	Indicates the manufacturer's batch code so that the batch or departure can be identified.
	Catalog number	Indicates the manufacturer's reference number so that the product for the health can be identified.
	Sterilized using irradiation	Indicates that the product is sterilized
	Non-Sterile	Indicates that the product is not sterile
	Recyclable	Indicates that packaging can be recycled
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	Authorized representative in the European Community	Indicates the Authorized representative in the European Community
	CE mark	Indicates that the product has the CE marking and the identification of the Notified Body.
	Keep away from sunlight	Indicates a medical device that requires protection against light sources.

This review cancels and replaces the previous revision.

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