Instruction for use

POLYGLYCOLIC ACID FAST

(BRAIDED AND COATED POLYGLYCOLIC ACID - FAST ABSORBABLE)
STERILE SYNTHETIC ABSORBABLE SURGICAL SUTURE

DESCRIPTION

Braided and Coated Polyglycolic Acid - Fast Absorbable suture (Polyglycolic acid), is a sterile synthetic multifilament braided fast absorbable surgical suture prepared from homopolymers of 100% Glycolide. Braided and Coated Polyglycolic Acid - Fast Absorbable suture uses a polymer material with a lower molecular weight than regular Braided and Coated Polyglycolic Acid (Polyglycolic acid suture), making it fast absorbable leading to rapid loss of tensile strength after implantation. Braided and Coated Polyglycolic Acid - Fast Absorbable suture is available in undyed form. For added lubrication, smooth passage through tissue and easy knot tie down, Braided and Coated Polyglycolic Acid - Fast Absorbable suture is coated with a mixture of polycaprolactone and calcium stearate. The composition of polymer material Polyglycolic acid used to prepare the suture is ≥ 95% and the composition of coating materials polycaprolactone and calcium stearate is ≤ 5%. Braided and Coated Polyglycolic Acid - Fast Absorbable suture has been found to be non-antigenic, non-pyrogenic and elicits only a mild to moderate tissue reaction during the process of absorption. Braided and Coated Polyglycolic Acid - Fast Absorbable sutures comply with the requirements established by the European Pharmacopoeia (E.P) for Sterile Synthetic Absorbable Braided Sutures, and United States Pharmacopeia (U.S.P) for Absorbable Surgical Suture (Synthetic).

INTENDED USE

Braided and Coated Polyglycolic Acid - Fast Absorbable suture is intended for use in general soft tissue approximation where only short term wound support is required and where rapid absorption of the suture would be beneficial. Due to the rapid absorption profile, Braided and Coated Polyglycolic Acid - Fast Absorbable suture is useful for skin closure, particularly in paediatric surgery, episiotomies, circumcision and closure of oral mucosa. Braided and Coated Polyglycolic Acid - Fast Absorbable suture is also useful for conjunctival suturing in ophthalmic surgery.

INDICATION

Indicated in surgical procedures where a combination of an absorbable suture and short term wound support (Upto 7 days) is required for wound closure and soft tissue approximation.

CONTRAINDICATIONS

Due to its rapid loss of tensile strength, Braided and Coated Polyglycolic Acid - Fast Absorbable suture should not be used where extended approximation of tissue under stress is required or where wound support beyond 7 days is required. Braided and Coated Polyglycolic Acid - Fast Absorbable suture should not be used for ligation, and in cardiovascular and neurological tissues.

INTENDED USERS

Trained and registered health care professional only.

INTENDED PATIENT POPULATION

This suture can be used in patients irrespective of age and sex, in line with the intended use, indications and contraindications.

PERFORMANCE

Braided and Coated Polyglycolic Acid - Fast Absorbable suture after implantation elicits a mild to moderate initial inflammatory reaction in tissues and is eventually replaced with an ingrowth of fibrous connective tissue. Braided and Coated Polyglycolic Acid - Fast Absorbable suture being absorbable in nature gets absorbed in the tissues by means of hydrolysis, during which the Polyglycolic acid polymer is degraded into Glycolic acid which is eventually absorbed and metabolized in the body. This process of absorption initially leads to progressive loss of tensile strength, followed by a loss of mass.

Implantation studies in rats show the following profile:

Days post Implantation	Approximate % Original Strength Remaining
7 days	44%
14 days	14.6%
21 days	All original strength loses within 21 days

Absorption of Braided and Coated Polyglycolic Acid - Fast Absorbable suture is essentially complete approximately between 42 to 63 days.

APPLICATION

Sutures should be selected and implanted depending on patient's condition, surgical experience, surgical technique and wound size. Braided

and Coated Polyglycolic Acid - Fast Absorbable typically falls off in seven to ten days post operatively or can be wiped off subsequently with sterile gauze. Normally the removal of the sutures is not required.

WARNINGS

- This suture is intended to be used by a trained and registered health care professional.
- Users should be familiar with surgical procedures and surgical suturing techniques involving absorbable sutures before using Braided and Coated Polyglycolic Acid - Fast Absorbable suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Users should consider the in vivo performance when selecting a suture.
- Users must apply their professional judgment when determining the appropriate sized suture based on the specific indication, wound size, preferred surgical technique, patient's condition and history.
- The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing.
- As with any foreign body, prolonged contact of this suture with salt solutions such as those found in the urinary or biliary tracts may result in calculus formation.
- As an absorbable suture, it may transiently act as a foreign body and as with any foreign body, presence of bacterial contamination may enhance bacterial infectivity. Thus, acceptable surgical practice should be followed with respect to drainage and closure of infected or contaminated wounds.
- Wound dehiscence may occur when suture fails to provide adequate wound support in closure of the sites where expansion, stretching or distension occur. Hence use of supplemental non-absorbable suture should be considered by the user in the closure of sites which may undergo expansion, stretching, or distension, or which may require additional support.
- Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the user.
- Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.
- Conjunctival, cuticular and vaginal epithelial sutures remaining in place for extended periods may be associated with localized irritation and should be snipped off or removed as indicated.
- Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.

PRECAUTIONS

- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance & cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage of the needle. Reshaping needles may cause them to lose strength and make less resistant to bending and breaking. Broken needles may result in extended or additional surgeries or residual foreign bodies.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard used needles in 'Sharps' containers.
- When handling this suture, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- When required, straighten sutures with a gentle pull, never crush or rub them.
- As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the user.
- Do not tie bands around pouches/packs of sutures as it may damage the pouches/pack.
- Do not re-sterilize, do not reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

SIDE EFFECTS

Side effects associated with the use of Braided and Coated Polyglycolic Acid - Fast Absorbable suture include allergic response in certain patients,

transient local irritation at the wound site, transient inflammatory foreign body response, mild to moderate inflammatory tissue reaction, erythema and induration during the absorption process of subcuticular sutures, and rarely wound dehiscence. Like all foreign bodies, Braided and Coated Polyglycolic Acid - Fast Absorbable suture may enhance an existing infection.

CLINICAL BENEFITS:

The clinical benefits of this suture are, its lubricant coating creates a smooth surface thereby decreasing the coefficient of friction allowing easy tissue passage and ensuring minimal/no tissue drag and trauma. It being braided provides secure knot holding capacity, when compared to monofilament sutures. When used for tissue approximation following injury or surgery, it helps in closing dead space, minimizes blood loss and infection, and also aid in wound healing. When used for approximating skin edges, it provides an aesthetically pleasing and functional result. It provides adequate initial tensile strength as per the U.S.P and E.P requirements. It provides adequate short term wound support during the critical wound healing period, even with smaller diameter sutures. It has a predictable absorption rate enabling decision making regarding the application of sutures in relevant tissues and conditions, as against natural sutures like catgut which doesn't have a predictable absorption rate. It is absorbed by the process of hydrolysis enabling moderate tissue reaction and inflammatory response. This suture being fast absorbable in nature, eliminates the discomfort of removal of sutures, unlike non-absorbable sutures thereby saving the surgeon's time and lessening patient anxiety and discomfort. This suture is mainly beneficial when used for suturing sensitive areas in procedures like episiotomy and circumcision, as they provide adequate short term wound support during the critical healing period, and then suture gets absorbed/falls off on their own without the need for suture removal. Also when used in episiotomy repair, it leads to decreased incidence of dyspareunia, which will positively improve the quality of life in patients.

STERILITY

Braided and Coated Polyglycolic Acid - Fast Absorbable suture is sterilized by ethylene oxide gas. Do not re-sterilize, do not reuse. Do not use if package is opened or damaged. Discard opened, unused sutures.

SUPPLY

Braided and Coated Polyglycolic Acid - Fast Absorbable suture is available in various U.S.P. sizes and E.P Metric sizes. number. Sutures are supplied sterile in pre-cut lengths, both non-needled and attached to various needle type, shape and length. Sutures are packed in a printed box with quantity indicated on the box label.

STORAGE

Recommended storage conditions: Below 30° C, away from moisture and direct sunlight. Do not use after expiry date.

DISPOSAL

Discard used sutures and needles contaminated with blood and tissue in the container meant for infectious waste. Unused expired pouches should be incinerated or disposal should be done as per local regulations.

SYMBOLS USED ON LABEL

Do not re-use

Upper limit of temperature

Batch code

Caution

Mfg. Date & Country of manufacture

Use-by date

Keep away from sunlight

Manufacturer

Do not re-sterilize

Medical Device

Keep dry

Consult instructions for use

Sterilized using Ethylene oxide

Single sterile barrier system with protective packaging inside

Double sterile barrier system

Do not use if package is damaged

Authorized representative in the European Union

Importer

Distributor

Unique device identifier

GS 1 Data matrix barcode

Catalogue number

European conformity mark and Identification of Notified Body. The product meets the safety and performance requirements of EU Medical devices Regulation (EU) 2017/745

BASIC UDI-DI : 8903837H0101010201ZH

Healthium Medtech Limited
#472-D, 13th Cross, 4th Phase,
Peenya Industrial Area, Bangalore, Karnataka-560058, India
Email: care@healthiummedtech.com
Customer care number: +91-80-41868198
Mfg. Lic. No.: MFG/MD/2021/000367

WAM
120 rue Jean de Guirmand
13290 Aix-en-Provence
France

MED DEVICES LIFESCIENCES B.V.
Keizersgracht 482,
1017 EG Amsterdam, The Netherlands
Email: info@meddevices.net
Phone: +31-202254558

Summary of safety and clinical performance(SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the BASIC UDI-DI : 8903837H0101010201ZH
<https://ec.europa.eu/tools/eudamed>

Adverse events should be reported to Healthium Medtech Ltd, or an in-country representative, and to the health authority where the event occurred. Email id: care@healthiummedtech.com

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