POLYPROPYLENE
(MONOFILAMENT POLYPROPYLENE)
STERILE NON-ABSORBABLE SURGICAL SUTURE

DESCRIPTION:

Monofilament Polypropylene suture (Polypropylene) is a sterile synthetic monofilament non-absorbable surgical suture prepared from a synthetic linear Polyolefin, Polypropylene. The molecular formula is (C₃H₆)_n. Monofilament Polypropylene suture is coloured blue, dyed with Phthalocyanine blue (conforming to U.S.Code of Federal regulations). Monofilament Polypropylene suture may also be supplied with PTFE (polytetrafluoroethylene) pledgets for use as a pad between the suture and the tissue surface to increase the load bearing area. The composition of Polypropylene suture is ≥ 99.5 % of Polypropylene and ≤ 0.5% of Phthalocyanine blue. Monofilament Polypropylene sutures comply with the requirements established by the European Pharmacopoeia (E.P) for Sterile Non-Absorbable sutures, and United States Pharmacopeia (U.S.P) for Non-Absorbable Surgical Suture.

INTENDED USE:

Monofilament Polypropylene suture is intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological tissues.

INDICATION:

Indicated for skin closure, and in surgical procedures where a combination of a non-absorbable suture and long term wound support is required for soft tissue approximation and/or ligation.

CONTRA INDICATION:

None known.

INTENDED USERS:

Trained and registered health care professional only.

INTENDED PATIENT POPULATION:

This suture can be used in patients irrespective of age and sex, in line with the intended use, indications and contraindications.

PERFORMANCE:

Monofilament Polypropylene suture being a monofilament suture is remarkably smooth and easy to handle. It does not have interstices in its structure and being a monofilament suture, has been successfully employed in surgical wounds which could subsequently become infected or contaminated where it can minimize later sinus formation and suture extrusion. It elicits a minimal initial inflammatory reaction in tissues which is followed by gradual encapsulation of the suture by fibrous connective tissue. It is neither absorbed, nor it is subjected to degradation or weakening by the action of tissue enzymes. Due to its relative biological inertness, it is recommended for use in tissues, where the least possible suture reaction is desired. Because of its lack of adherence to tissues, it is effective as a pull out suture.

APPLICATION:

Suture should be selected and implanted depending on patient's condition, surgical experience, surgical technique and wound size.

WARNINGS:

- This suture is intended to be used by a trained and registered health care professional.
- Users should be familiar with surgical procedures and surgical suturing techniques involving non-absorbable sutures before using Monofilament Polypropylene suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.
- Users must apply their professional judgment when determining the appropriate sized suture based on the specific indication, wound size, preferred surgical technique, patient's condition and history.
- As with any foreign body, prolonged contact of this suture with salt solutions such as those found in the urinary or biliary tracts may result in calculus formation.
- As a non-absorbable suture, it may act as a foreign body and as with any

foreign body, presence of bacterial contamination may enhance bacterial infectivity. Thus, acceptable surgical practice should be followed with respect to drainage and closure of infected or contaminated wounds.

- Wound dehiscence may occur when suture fails to provide adequate wound support in closure of the sites where expansion, stretching or distension occur. Postoperatively and until healing is complete, fixation provided by this suture should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this suture should be protected. The postoperative regimen prescribed by the surgeon should be strictly followed to avoid adverse stresses applied to the suture.
- Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the user.
- Device removal when required, should be followed by adequate postoperative management.

PRECAUTIONS:

- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance & cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage of the needle. Reshaping needles may cause them to lose strength and make less resistant to bending and breaking. Broken needles may result in extended or additional surgeries or residual foreign bodies.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard used needles in 'Sharps' containers.
- When handling this suture, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- When required, straighten sutures with a gentle pull, never crush or rub them.
- As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the user. The use of additional throws may be particularly appropriate when knotting a monofilament suture.
- Do not tie bands around pouches/packs of sutures as it may damage the pouches/pack.
- Do not re-sterilize, do not reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

SIDE EFFECTS:

Side effects associated with the use of Monofilament Polypropylene suture include allergic response in certain patients, transient local irritation at the wound site, minimal initial inflammatory tissue reaction, pain, edema, and erythema at the wound site, wound dehiscence, and calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine or bile occurs. Like all foreign bodies, Monofilament Polypropylene suture may enhance an existing infection.

CLINICAL BENEFITS:

The clinical benefits of this suture are, it being a monofilament suture allows easy tissue passage and ensures minimal/no tissue drag and trauma, and because of its lack of adherence to tissues, it is effective as a pull out suture. It provides secure knot holding capacity, and being monofilament in nature doesn't have interstices in its structure leading to reduced incidence of bacterial colonization, when compared to braided non-absorbable sutures. When used for tissue approximation and/or ligation following injury or surgery, it helps in closing dead space, minimizes blood loss and infection, and also aid in wound healing. It provides adequate tensile strength as per the U.S.P and E.P requirements. It provides adequate wound support during the critical wound healing period, even with smaller diameter sutures. Owing to the availability of very small sizes and with adequate tensile strength, it will be very useful and relevant in ophthalmic and cosmetic surgery also. It exhibits all the desirable physical characteristics of the ideal suture material for vascular anastomosis like long-term durability, high

tensile strength, a favourable stress-strain relationship, minimal biological reactivity and infection, a low coefficient of friction, and knot security making it an ideal suture material for artery anastomosis. It provides a long term tensile strength retention enabling decision making regarding the application of sutures in relevant tissues and conditions, as against absorbable sutures which provides only short term wound support. It is biocompatible with mild tissue reaction and inflammatory response. Due to its relative biological inertness, it is recommended for use in tissues, where the least possible suture reaction is desired. The pain associated with removal of sutures is lesser when compared to braided non-absorbable sutures and skin staples. This suture leads to less reactive scars with a lower tendency towards hypertrophic scar formation, and hence provides good cosmetic results, and therefore very beneficial for closing skin incisions to provide a superior cosmetic outcome. It has long term durability and it spreads tension evenly across the wound edges without causing skin edge necrosis and most suitable for delayed primary closure of fasciotomy wounds and mass closure of vertical laparotomy wounds.

STERILITY:

Monofilament Polypropylene suture is sterilized by ethylene oxide gas. Do not re-sterilize, do not reuse. Do not use if package is opened or damaged. Discard opened, unused sutures.

SUPPLY:

Monofilament Polypropylene suture is available in various U.S.P. sizes and E.P Metric sizes. Sutures are supplied sterile in pre-cut lengths, both non-needed and attached to various needle type, shape and length. Sutures are packed in a printed box with quantity indicated on the box label.

STORAGE:

Recommended storage conditions: Below 30° C, away from moisture and direct sunlight. Do not use after expiry date

DISPOSAL:

Discard used sutures and needles contaminated with blood and tissue in the container meant for infectious waste. Unused expired pouches should be incinerated or disposal should be done as per local regulations.

SYMBOLS USED ON LABEL:

Do not re-use

30° C Upper limit of temperature

Batch code

Caution

Mfg. Date & Country of manufacture

Use-by date

Keep away from sunlight

Manufacturer

Do not re-sterilize

Medical Device

Keep dry

Consult instructions for use

Sterilized using Ethylene oxide

Single sterile barrier system with protective packaging inside

Do not use if package is damaged

Authorized representative in the European Union

Importer

Distributor

Unique device identifier

GS 1 Data matrix barcode

Catalogue number

European conformity mark and Identification of Notified Body. The product meets the safety and performance requirements of EU Medical devices Regulation (EU) 2017/745

BASIC UDI-DI : 8903837H010201010426

Healthium Medtech Limited
#472-D, 13th Cross, 4th Phase,
Peenya Industrial Area, Bangalore, Karnataka-560058, India
Email: care@healthiummedtech.com
Customer care number: +91-80-41868198
Mfg. Lic. No: MFG/MD/2021/000367

WAM

WAM
120 rue Jean de Guiramand
13290 Aix-en-Provence
France

MED DEVICES LIFESCIENCES B.V.

Keizersgracht 482,
1017 EG Amsterdam, The Netherlands
Email: info@meddevices.net
Phone: +31-202254558

Summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the BASIC UDI-DI : 8903837H010201010426
<https://ec.europa.eu/tools/eudamed>

Adverse events should be reported to Healthium Medtech Ltd, or an in-country representative, and to the health authority where the event occurred. Email id: care@healthiummedtech.com

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