Instruction for use

SILK
(BLACK BRAIDED SILK)
STERILE NON-ABSORBABLE SURGICAL SUTURE

DESCRIPTION:

SILK suture is a sterile multifilament braided non-absorbable surgical suture prepared from an organic protein called fibroin. This protein is derived from the domesticated species Bombyx mori of the family Bombycidae. Natural Silk contains approximately 30% extraneous material or gum which is removed by "Degumming" process. This process is essential for a compact braid while ensuring that the filaments retain their natural body and elasticity. SILK suture is dyed with black Hematoxylin and coated with bees wax. The composition of Black braided silk suture is 84%-94% of Silk fibroin, 5%-15% of Bees wax, and ≤1% of Hematoxylin dye. SILK sutures comply with the requirements established by the European Pharmacopoeia (E.P) for Sterile Non-Absorbable sutures, and United States Pharmacopoeia (U.S.P) for Non-Absorbable Surgical Suture.

INTENDED USE:

SILK suture is intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological tissues.

INDICATION:

Indicated for skin closure, and in surgical procedures where a combination of a non-absorbable suture and long term wound support (except where permanent retention of tensile strength is required) is required for soft tissue approximation and/or ligation.

CONTRA INDICATION:

Use of this suture is contra indicated in patients with known sensitivities or allergies to Silk. Due to gradual loss of tensile strength which may occur over prolonged periods in vivo, SILK suture should not be used where permanent retention of tensile strength is required. For example: as in fixation of intraocular lenses or synthetic vascular grafts. Do not expose the suture pack to chemical disinfectants containing oxidizing agents like Hydrogen peroxide or other similar chemicals, as such exposure may affect the product quality.

INTENDED USERS:

Trained and registered health care professional only.

INTENDED PATIENT POPULATION:

This suture can be used in patients irrespective of age and sex, in line with the intended use, indications and contraindications.

PERFORMANCE:

SILK suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. While SILK suture is not absorbed, progressive degradation of the proteinaceous silk fibres in vivo may result in gradual loss of the suture's tensile strength over a period of time.

APPLICATION:

Suture should be selected and implanted depending on patient's condition, surgical experience, surgical technique and wound size.

WARNINGS:

- This suture is intended to be used by a trained and registered health care professional.
- Users should be familiar with surgical procedures and surgical suturing techniques involving non-absorbable sutures before using SILK suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.
- Users must apply their professional judgment when determining the appropriate sized suture based on the specific indication, wound size, preferred surgical technique, patient's condition and history.
- As with any foreign body, prolonged contact of this suture with salt solutions such as those found in the urinary or biliary tracts may result in calculus formation.
- As a non-absorbable suture, it may act as a foreign body and as with any foreign body, presence of bacterial contamination may enhance

bacterial infectivity. Thus, acceptable surgical practice should be followed with respect to drainage and closure of infected or contaminated wounds.

- Wound dehiscence may occur when suture fails to provide adequate wound support in closure of the sites where expansion, stretching or distension occur. Postoperatively and until healing is complete, fixation provided by this suture should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this suture should be protected. The postoperative regimen prescribed by the surgeon should be strictly followed to avoid adverse stresses applied to the suture.
- Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the user.
- Device removal when required, should be followed by adequate postoperative management.

PRECAUTIONS:

- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance & cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage of the needle. Reshaping needles may cause them to loose strength and make less resistant to bending and breaking. Broken needles may result in extended or additional surgeries or residual foreign bodies.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard used needles in 'Sharps' containers.
- When handling this suture, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- When required, straighten sutures with a gentle pull, never crush or rub them.
- As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the user.
- Do not tie bands around pouches/packs of sutures as it may damage the pouches/pack.
- Do not re-sterilize, do not reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

SIDE EFFECTS:

Side effects associated with the use of SILK suture include allergic response in certain patients who are known to be sensitive to silk, transient local irritation at the wound site, initial inflammatory tissue reaction, pain, edema, and erythema at the wound site, gradual loss of tensile strength over time, wound dehiscence, and calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine or bile occurs. Like all foreign bodies, SILK suture may enhance an existing infection.

CLINICAL BENEFITS:

The clinical benefits of this suture are, its lubricant coating creates a smooth surface thereby decreasing the coefficient of friction allowing easy tissue passage and ensuring minimal/no tissue drag and trauma. It being braided provides secure knot holding capacity, when compared to monofilament non-absorbable sutures. It is soft and very easy to handle. When used for tissue approximation and/or ligation following injury or surgery, it helps in closing dead space, minimizes blood loss and infection, and also aid in wound healing. When used for approximating skin edges, it provides a good functional result. It is relatively less expensive when compared with other currently available non-absorbable suture materials. It provides adequate tensile strength as per the U.S.P and E.P requirements. It provides adequate wound support during the critical wound healing period, even with smaller diameter sutures. It provides long term tensile strength retention enabling decision making regarding the application of sutures in relevant tissues and conditions, as against absorbable sutures which provides only short term wound support. Owing to this suture's minimum propensity to tear through

tissue it is very useful in gastrointestinal, dental, ocular and neural procedures. It is biocompatible with mild tissue reaction and inflammatory response. The pain associated with removal of sutures is lesser when compared to skin staples.

STERILITY:

SILK suture is sterilized by ethylene oxide gas. Do not re-sterilize, do not reuse. Do not use if package is opened or damaged. Discard opened, unused sutures.

SUPPLY:

SILK suture is available in various U.S.P. sizes and E.P Metric sizes. Sutures are supplied sterile in pre-cut lengths, both non-needed and attached to various needle type, shape and length. Sutures are packed in a printed box with quantity indicated on the box label. Non-needed SILK suture is also available as Non-sterile reels.

STORAGE:

Recommended storage conditions: Below 30° C, away from moisture and direct sunlight. Do not use after expiry date

DISPOSAL:

Discard used sutures and needles contaminated with blood and tissue in the container meant for infectious waste. Unused expired pouches should be incinerated or disposal should be done as per local regulations.

SYMBOLS USED ON LABEL:

	Do not re-use		30° C	Upper limit of temperature
	Batch code			Caution
	Mfg. Date & Country of manufacture			
	Use-by date			Keep away from sunlight
	Manufacturer			Do not re-sterilize
	Medical Device			Keep dry
	Consult instructions for use			
	Sterilized using Ethylene oxide			
	Single sterile barrier system with protective packaging inside			
	Do not use if package is damaged			
	Authorized representative in the European Union			
	Importer			Distributor
	Unique device identifier			
	GS 1 Data matrix barcode			
	Catalogue number			

 **2265** European conformity mark and Identification of Notified Body. The product meets the safety and performance requirements of EU Medical devices Regulation (EU) 2017/745

BASIC UDI-DI : 8903837H01020202032L

Healthium Medtech Limited
#472-D, 13th Cross, 4th Phase,
Peenya Industrial Area, Bangalore, Karnataka-560058, India
Email: care@healthiummedtech.com
Customer care number: +91-80-41868198
Mfg. Lic. No.: MFG/MD/2021/000367

WAM
120 rue Jean de Guiramand
13290 Aix-en-Provence
France

MED DEVICES LIFESCIENCES B.V.
Keizersgracht 482,
1017 EG Amsterdam, The Netherlands
Email: info@meddevices.net
Phone: +31-202254558

Summary of safety and clinical performance(SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the BASIC UDI-DI : 8903837H01020202032L
<https://ec.europa.eu/tools/eudamed>

Adverse events should be reported to Healthium Medtech Ltd, or an in-country representative, and to the health authority where the event occurred. Email id: care@healthiummedtech.com

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