正面

Instructions for use

【Product name】Mouth mirror
【Model】RS-22、RS-42、RS-52、RD-52、RX-Y03、RX-Y05、RX-F27、RX-F39、RX-Y10、RX-KIT、RV-SM02、RS-L01、RS-M-21、RS-M-41、RS-M-51、RD-M-51、RD-M-41、RV-FFM-01、RV-JH-2S、RV-JH-3S、RV-JH-4S、RV-JH-5S、RV-JH-4D、RV-SH-4S、RV-SH-5S、RV-SH-4D、RV-DP-01、RV-MM-G01、RV-MM-G146、RV-MM-G03
【Product description】Consists of a handle, a mirror with or without a connecting rod.Mirrors are made of stainless steel or polished glass. Not supplied sterile. Sterilized before use by the using body according to the instructions.
【Intended Use】For Oral examination.
【Applicable People】Adult and pediatric dental patients
【Method of use】
1、Please autoclave before first use and after each use.
2、Put the mouth mirror into the mouth, make sure the mirror surface is toward the teeth and gums, manipulate the handle with your fingers to move the mouth mirror in the mouth in order to check each area.
【Cleaning and Sterilization】Scrub the mirror with neutral detergent or hot water, then wrap and place the mouth mirror in an autoclave container for sterilization. Temperature (121℃±2℃) for 15 minutes or temperature (135℃±2℃) for 3.5 minutes.
【Cautions】
1、This product is for professional use only.
2、Check the mirror for cracks or other defects before use and replace it promptly if necessary.
3、When using the mouth mirror, keep your hand steady so that the mouth mirror does not slip off.
4、Don't insert the mouth mirror deep into the throat, so as not to cause breathing difficulties.
5、After using the mouthpiece, clean and sterilize it in time and store it to avoid bacterial growth.
【Contraindications】N/A
【Production date】See label or package.
【Storage Method】Keep at room temperature and keep away from light.
【Expiry date】 5 years.

背面

Labels, Packing Logo Design:

Symbol	Introductions	Symbol	Introductions
CE	CE Symbol	MD	Medical Device
	Name and Address of Manufacturer	EU REP	Name and Address of European Union Authorized Representative
	Date of manufacture	#	Model Number
CC	Country of Manufacture		Use-by Date
	Importer		Distributor
	Caution	UDI	Unique Device Identifier
LOT	Batch Code		Consult Instructions for Use
	Upper limit of temperaturelimit	134°C	Sterilizable in a steam sterilizer (autoclave) at temperature specified
	Keep Away From Sunlight		Keep Dry
NON STERILE	Non-sterile	REF	Catalogue number

Manufacturer Information

Company: Foshan Ruibang Xinsheng Medical Instrument Co., Ltd
Address: 501, 5/F, Building 9, Area D, Langsha New Light Industrial Base of Guangdong, Shishan Town,Nanhai District, Foshan City, Guangdong, 528220

European Authorized Representative

EU REP
SUNGO Europe B.V.
Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands

Issued date:2025-07-17

Ver: A/1